

# **UCLA**

## **UCLA Previously Published Works**

### **Title**

Drosophila hematopoiesis: Markers and methods for molecular genetic analysis

### **Permalink**

<https://escholarship.org/uc/item/42x9s8k2>

### **Journal**

Methods, 68(1)

### **ISSN**

1046-2023

### **Authors**

Evans, Cory J

Liu, Ting

Banerjee, Utpal

### **Publication Date**

2014-06-01

### **DOI**

10.1016/j.ymeth.2014.02.038

Peer reviewed

Published in final edited form as:

*Methods*. 2014 June 15; 68(1): 242–251. doi:10.1016/j.ymeth.2014.02.038.

## ***Drosophila* hematopoiesis: markers and methods for molecular genetic analysis**

**Cory J. Evans<sup>\*</sup>, Ting Liu<sup>\*</sup>, and Utpal Banerjee<sup>\*,§,\*\*,§§</sup>**

<sup>\*</sup>Department of Molecular, Cell and Developmental Biology, University of California, Los Angeles Los Angeles, CA 90095 USA

<sup>§</sup>Department of Biological Chemistry, University of California, Los Angeles Los Angeles, CA 90095 USA

<sup>\*\*</sup>Molecular Biology Institute, University of California, Los Angeles Los Angeles, CA 90095 USA

<sup>§§</sup>Eli and Edythe Broad Center of Regenerative Medicine and Stem Cell Research, University of California, Los Angeles Los Angeles, CA 90095 USA

### **Abstract**

Analyses of the *Drosophila* hematopoietic system are becoming more and more prevalent as developmental and functional parallels with vertebrate blood cells become more evident. Investigative work on the fly blood system has, out of necessity, led to the identification of new molecular markers for blood cell types and lineages and to the refinement of useful molecular genetic tools and analytical methods. This review briefly describes the *Drosophila* hematopoietic system at different developmental stages, summarizes the major useful cell markers and tools for each stage, and provides basic protocols for practical analysis of circulating blood cells and of the lymph gland, the larval hematopoietic organ.

### **Introduction**

Blood cells have been studied in invertebrate animal models for more than one hundred years, although how these cells relate to vertebrate blood cells in function and development has only recently begun at the molecular genetic level. This work is highlighted by studies in the fruit fly *Drosophila melanogaster*, an established premier model system for genetic studies in other contexts. The earliest studies of blood cells (known as hemocytes) in *Drosophila* date back to the late 1800s, however significant progress in this area did not begin in earnest until the 1950s when Rizki and Rizki began their studies of blood cell types and associated functions, particularly in the areas of innate immune responses and self-/non-self-recognition, and established the nomenclature for the *Drosophila* blood cell system that

© 2014 Elsevier Inc. All rights reserved.

**Corresponding author:** Utpal Banerjee Department of Molecular, Cell, and Developmental Biology University of California, Los Angeles Los Angeles, California, 90095, USA Phone: 1-310-206-5439 Fax: 1-310-206-3987 banerjee@mbi.ucla.edu.

**Publisher's Disclaimer:** This is a PDF file of an unedited manuscript that has been accepted for publication. As a service to our customers we are providing this early version of the manuscript. The manuscript will undergo copyediting, typesetting, and review of the resulting proof before it is published in its final citable form. Please note that during the production process errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

is still in use today [1, 2]. Furthermore, their observational work led to the earliest suggestion that the various *Drosophila* blood cell types arise from a common precursor cell type, an idea similar to that which would be established in the mouse within the next few years [3].

Since that time, with the advent of molecular genetics, quite a bit has been discovered about how *Drosophila* blood cells are specified and develop and how this relates to blood-forming processes in vertebrate systems. Much of this knowledge has been acquired with the identification and use of many new cellular markers, genetic tools, and analytical methods. This review will briefly describe the key elements of the *Drosophila* hematopoietic system and highlight useful reagents and approaches that are currently available.

## The *Drosophila* hematopoietic system

Mature *Drosophila* blood cells are found as at least three distinct types known as plasmatocytes, crystal cells and lamellocytes. The plasmatocyte behaves as a macrophage-like cell and is the predominant cell type comprising more than 90% of the hemocyte repertoire [1, 4]. Crystal cells, named for their cytoplasmic paracrystalline protein inclusions, represent approximately 5% of hemocytes in circulation while lamellocytes are rarely observed. Lamellocytes can be induced to differentiate, however, in the context of various conditions including immune challenge by parasitic wasps [5-7]. Several additional cell types have been described based on morphological features or as precursor populations lacking differentiation markers [2, 8-10], but their identification as independent cell types will require the discovery of new and unique markers.

A major function of hemocytes is in the provision of cellular innate immune responses, which is achieved primarily through the phagocytic clearance of microbial pathogens or the encapsulation of larger parasites. Additionally, plasmatocytes are known to promote the humoral immune response by secreting cytokine-like proteins and antimicrobial peptides [11-15]. Accordingly, flies with impaired plasmatocyte function are more susceptible than wild-type flies to microbial infection and are less effective in the potentiation of humoral immune responses [16-19]. Hemocytes also have roles during development where they secrete and remodel extracellular matrix components as well as remove cellular debris by phagocytosis [20-23]. These roles are highlighted by their requirement in the embryonic nervous system, the development of which is severely disrupted in the absence of plasmatocyte function [24].

As professional phagocytes, plasmatocytes use a large repertoire of receptors to identify and engulf different targets. Receptors recognizing dying cells include Croquemort (Crq), a CD36 family member [25], PSR, a phosphatidyl serine receptor [26], and Draper, which is homologous to the CED-1 phagocytosis receptor in *C. elegans* [27, 28]. Numerous receptors that bind microbial pathogens are also known and include the peptidoglycan recognition protein PGRP-LC and the scavenger receptor D-SR-CI [29-31]. Additionally, the receptors Eater, Nimrod (NimC1, also known as the P1 antigen) and Dscam have been found to be important in bacterial clearance. Eater and Nimrod both contain multiple EGF-like motifs, called NIM repeats, which, in the case of Eater, have been shown to bind bacteria and

mediate an array of intermolecular interactions [32, 33]. Nimrod belongs to a family of ten related proteins (grouped into three classes, NimA-NimC), the majority of which are expressed by hemocytes [33]. Interestingly, the NimB proteins appear to function as opsonins because they enter the secretory pathway but lack transmembrane and cytosolic domains. Several isoforms of Dscam have also been identified in the hemolymph suggesting that it may also function to opsonize microbes [34]. The extracellular domain of Dscam receptors is composed of repeating immunoglobulin-like domains, a subset of which are highly variable due to alternative splicing [35]. Though microarray analysis has shown that more than 18,000 alternatively-spliced isoforms are expressed by hemocytes, just five isoforms comprise 80-90% of Dscam mRNAs in these cells, suggesting a blood-specific role for these particular Dscam types [34].

While plasmatocytes primarily function in phagocytic clearance, the main function of crystal cells is in the process of melanization, the darkening and hardening of tissue due to the local deposition of melanin. Melanization commonly occurs during immune responses in the context of barrier formation around pathogens too large to remove by phagocytosis, as well as during the process of wound healing, particularly when the cuticular epithelia has been breached. Phenoloxidase (PO) enzymes mediate the oxidation of phenols into quinones that then polymerize into melanin. Tightly-regulated serine protease cascades convert Prophenoloxidase (PPO) zymogens into the active PO form [36-39]. Experimental evidence indicates that the large, cytoplasmic inclusions that crystal cells exhibit (and from which they derive their name) are composed primarily of PPO [2, 5, 10, 40]. Furthermore, *Drosophila* mutant lines that lack crystal cells also lack PO activity in the circulating hemolymph, identifying crystal cells as the primary source of this activity [40-43]. Genetic analysis suggests that crystal cell rupture is the mechanism by which PO activity is delivered to the hemolymph, a process mediated by JNK signaling and the TNF homolog Eiger [41]. Despite the distribution of PPO throughout the hemolymph, PO activity is tightly regulated and spatially restricted, such as at sites of coagulation and clot formation [44].

Lamellocytes are the most morphologically distinct blood cell type, being large (15-40  $\mu\text{m}$  across), disc-shaped cells. These cells function during the encapsulation response, where a cellular barrier, which includes plasmatocytes and crystal cells [45], forms around foreign objects that cannot be removed by phagocytosis. Normally, very few lamellocytes are observed [46]; however numerous lamellocytes can be induced to differentiate in response to signals that include parasitic wasp infestation, injection of foreign objects into the hemocoel, and sterile wounding [5, 6, 13, 47].

With regard to understanding innate immune mechanisms, *Drosophila* has proven to be an exceptional model that has provided key insights into systems, such as Toll-like receptor signaling, with relevance to human biology [48]. The majority of this work in flies has focused on understanding the humoral response to microbial pathogens, however the role of hemocytes in the immune context is beginning to be explored in more detail. Methodologies and protocols describing microbial challenge and the monitoring immune signaling have been previously published by Lemaitre and colleagues [49], (Editor: please add current reference if Lemaitre contributes to this issue). A detailed protocol demonstrating wasp parasitization and associated analytical methods has also recently been made available [50].

Much like in vertebrate systems, hematopoietic development in *Drosophila* is biphasic with regard to timing and location [8, 51, 52]. *Drosophila* blood cells are first specified in the head mesoderm and migrate throughout the developing embryo [53]. At later stages of embryogenesis, a second group of hematopoietic precursors are specified within the cardiogenic mesoderm (thoracic segments 1-3) and form the lymph gland, a specialized organ that supports the hematopoietic process throughout subsequent larval development [54, 55]. At the end of the larval stage, the lymph gland breaks down and releases blood cells into circulation [56, 57]. Several reports have also suggested that hematopoietic cell populations exist as sessile cells in association with the internal larval body wall and/or the dorsal vessel (heart), however a better understanding of the blood-forming potential of these cells will require further exploration [58-60]. It has been demonstrated that circulating blood cells derived from both the embryonic head mesoderm and the lymph gland persist into the adult stage [56], however a major unanswered question is whether any adult blood cells arise from *de novo* hematopoiesis during the pupal or adult stages.

## Analysis of embryonic hemocytes

Molecular genetic analysis of embryonic blood cells has made use of standard *in situ* hybridization and immunostaining methods in combination with specific markers. The number of genes with expression patterns overlapping with embryonic hemocytes has grown in the last several years, however with regard to early developmental analysis of blood cells derived from the head mesoderm, several markers stand out: serpent (*srp*), glial cells missing (*gcm*), lozenge (*lz*), Prophenoloxidase A1 (*ProPOA1*), collagen (*Cg25c*), and Peroxidase (*Pxn*; see Table 1). The onset of *Srp* (a GATA transcriptional regulator) expression in the head mesoderm defines blood cell identity, while the expression of *Gcm* and *Lz* (a Runx family transcriptional regulator) mark the plasmatocyte and crystal cell lineages, respectively [8, 52]. Maturing macrophages express *collagen* and *Pxn*, while *ProPOA1* is a marker of mature crystal cells. The lamellocyte cell type has not been reported in the embryo. Each of the described markers has been analyzed by *in situ* hybridization and/or antibody staining, and some excellent examples describing their specific use, either alone or in combination, can be found in Waltzer et al., (2003) [61] and Milchanowski et al., (2004) [62], among others.

Embryonic lymph gland cells express *srp*, similar to head mesoderm blood cells, however the lymph gland lineage also expresses *odd skipped* (*odd*) [54, 63]. Because of the existence of *Odd* antibody and the reporters *odd-lacZ* and *odd-gal4*, *odd* expression has been a marker of choice for identifying lymph gland cells in late-stage embryos. Expression of the homeotic gene *Antennapedia* (*Antp*) also defines the earliest subdivision (stage 11) of the lymph gland by specifying Posterior Signaling Center (PSC) cells, which behave as niche-like cells at later developmental stages [64]. The EBF factor Collier (*Col*, also known as *Knot*) is also an important regulator and marker of the PSC (stage 11 onward) [63]. Collier functions downstream of *Antp* in the PSC and is important for its maintenance during larval development [64]. In *col* mutant animals, the PSC is lost, which causes lymph gland progenitors to differentiate prematurely [63]. Additional embryonic lymph gland markers are described in Table 1.

*Drosophila* embryos have also been ideal for research groups interested in live imaging of blood cells. Because of their small size, relative transparency, and the vast genetic tools that are available, live imaging of embryonic hemocytes has been used to address several basic developmental and cell biological questions, such as cell migration mechanics and chemotaxis [65-68]. Because of the versatility of this system, relatively recent detailed protocols describing live imaging and *in vivo* tracking of embryonic hemocytes have been published [69, 70].

## Analysis of larval hemocytes

Over the last several years, analysis of *Drosophila* hematopoiesis has primarily focused on the larval stages. Recently, there has been growing interest in characterizing larval circulating cells in both functional and hematopoietic contexts, and several protocols for *in vitro*, *ex vivo*, and *in vivo* hemocyte analysis have been described [58, 70-76]. A basic protocol for the collection and preparation of circulating larval hemocytes for analysis is described in Box 1. Most larval analyses, however, have examined blood development in the lymph gland, which has proven to be a useful system for understanding mechanisms of progenitor cell maintenance and differentiation [8, 77]. Although the lymph gland is specified in the embryo, its growth and differentiation into mature blood cells occurs during the larval stages. Owing to a lack of good molecular markers, early work was limited mainly to observation. However, many new markers and tools have been identified (see Table 2) that, in combination with the application of advanced microscopic techniques, have greatly enhanced the current understanding of how hematopoiesis occurs in the lymph gland. In particular, these advances demonstrated that the primary lobes of the lymph gland can be subdivided into distinct cellular populations that are spatially organized [54]. The periphery of each primary lobe contains maturing blood cells was named the Cortical Zone (CZ), whereas the juxtaposed medial region contains blood progenitor cells and was termed the Medullary Zone (MZ). The posterior tip of each primary lobe harbors the Posterior Signaling Center (PSC), several dozen specialized lymph gland cells that do not behave as blood cells but rather form a niche-like group that supports hematopoietic development [64, 78, 79].

Among lymph gland markers there are several that stand out as being particularly useful for genetic analysis (see Table 2 for these and others). The first lymph gland zone to be discovered was the PSC through the expression of Serrate (*Ser9.5-lacZ*), encoding a Notch receptor ligand. The *Serrate* enhancer, which is active in the PSC from the second instar onward, was subsequently utilized to generate a Gal4 driver line utilizing a fluorescent reporter (*UAS-GFP*) and as a genetic tool to manipulate PSC cells [54, 79]. Later research demonstrated that *Antennapedia*, as described previously, is expressed in and specifies PSC cells [64]. *Antennapedia* expression in the PSC is maintained throughout lymph gland development, making it a marker and driver of choice for PSC-related experiments (both *Antennapedia-gal4* and a monoclonal antibody are available, see Table 2). Similarly useful reagents are antibodies against Collier as well as the *col-gal4* reporter line, both of which mark PSC cells throughout larval stages [78, 80]. PSC cells also express and secrete Hedgehog protein, which is sensed by medullary zone progenitor cells and is important for their long-term maintenance within the lymph gland [64]. A direct *hedgehog-GFP* reporter

is available that allows for direct visualization of PSC cells [81], and can be placed in the context of other Gal4 drivers for analysis of non-autonomous PSC regulation (something not possible when using *Antennapedia-gal4* as a reporter, for example).

The medullary and cortical zones were identified simultaneously through their mutually exclusive expression of progenitor and mature cell markers, respectively. The classic medullary zone marker is *domeless-gal4* (*domeless* encodes a receptor that activates the JAK/STAT pathway), which was first isolated as an enhancer trap [54, 82]. Subsequently the *domeless* enhancer responsible for medullary zone expression was identified and has been used as a Gal4-independent *lacZ* (*domeless-MESO-lacZ*) reporter [78, 83]. Recent analysis demonstrating lymph gland regulation by the central nervous system (CNS) found that *domeless-gal4* is also expressed by various neurons in the brain. To circumvent CNS effects when using *domeless-gal4* to manipulate medullary zone blood progenitors, *ELAV-gal80* has been placed in combination with *domeless-gal4* [84]. The *ELAV* enhancer is neuron-specific and expresses high levels of Gal80, a natural protein inhibitor of Gal4, thereby mitigating any brain effects due to *domeless-gal4* expression. Additional useful genetic markers for lymph gland medullary zone progenitors include E-cadherin [54], Cubitus interruptus [64], Wingless [85], *TepIV-gal4* [86], phospho-Akt (p-Akt) [87], and Bag of marbles [88] (see Table 2).

Previous analysis of mitochondrial function in the lymph gland demonstrated that blood progenitor cells exhibit relatively high levels of reactive oxygen species (ROS), and that these ROS establish a critical signaling threshold for proper progenitor maintenance and differentiation [89]. Since then, fluorescent staining for ROS levels has become a standard analytic procedure for assessing how genetic change or immune challenge affects ROS within lymph glands and circulating cells. Many oxidation-sensitive dyes are available, however dihydroethidium (DHE) has been the marker of choice because of its high specificity for superoxide radicals and its ability to freely permeate cell membranes. Upon oxidation, fluorescent DHE metabolites are well retained by cells and tolerate mild fixation [90], which may be critical for simultaneous analysis of other markers such as GFP. Brief protocols for DHE analysis in both lymph glands and circulating blood cells are included here (see Boxes 5 and 6), with a more detailed protocol available as an online resource [91].

Differentiating cells of the lymph gland cortical zone express several useful markers, most notably *Hemolectin-gal4* and the P1 antigen, which are both blood specific [54]. The *Hemolectin-gal4* driver is one of the earliest known markers defining the onset of progenitor differentiation within the lymph gland at the mid-second instar, and stays on in mature cells [54, 92]. Expression of Peroxidase (both the protein and the *Peroxidase-gal4* reporter) [93, 94] is also an early marker of cortical zone formation [54], however it is also expressed at various levels in other tissues such as the fat body and brain. The *Collagen-gal4* driver [95] is slightly later than Peroxidase expression in the cortical zone [54], but also has strong expression in the fat body and, at lower levels, in various cell types [95]. The P1 antigen is specific to mature plasmotocytes [33, 96] and, by comparison, is considered to be a relatively late marker [54]. The P1 antigen was identified as one of several proteins interacting with monoclonal antibodies derived from *Drosophila* blood cell preparations [96]. The P1 antigen was subsequently identified as the protein product of the *Nimrod C1*

gene (*NimC1*), described previously [33]. While the anti-P1 monoclonal antibody remains among the first in the choice of tools to be used for blood analysis given its specific expression in late plasmatocytes, ease of use, and widespread availability, it should be noted that several commonly used *Drosophila* parental stocks are P1-negative [97]. This recessive expression polymorphism can lead to serious problems with experimental interpretation of progeny phenotypes, making it imperative to analyze P1-expression in parental strains. The P1-negative genetic background of any important experimental stock can easily be corrected through standard chromosome mechanics [97].

Crystal cells in the cortical zone can be identified by expression of the determinant gene *lozenge* (both the protein and the *lozenge-gal4* reporter) and mature marker Prophenoloxidase A1 (ProPOA1; protein). The gene encoding ProPOA1 is the *Black cells* (*Bc*) gene, and the crystal cell enhancer has been used to generate several different fluorescent, *Gal4*-independent *Bc*-expression reporters [81]. Lamellocytes also now have several useful markers available. A classic marker of lamellocyte differentiation is a misshapen (a JNK activator)-*lacZ* reporter gene (*msn-lacZ*) [98], initially isolated as an enhancer trap. The *msn* lamellocyte enhancer has subsequently been identified and used to generate several useful *Gal4*-independent fluorescent reporter lines [99]. Useful lamellocyte protein markers, for which antibodies are available, include L1 (identified along with P1 described above) [96], Filamin-240 [100],  $\alpha$ -PS4 integrin [78], and Myospheroid  $\beta$ PS-integrin [101]. A GFP enhancer trap of the L1 gene (*atilla*) is also available [102].

## Development of new tools for lymph gland analysis

Although several new markers for lymph gland analysis have been identified, their relative utility varies. As described, some markers are antibody based (many of which are not monoclonal) and their availability and quality can vary widely. Genetic reporters such as *Hml-gal4*, while extremely useful, can be difficult for experiments in which it is necessary to discern whether *Gal4* is a cell identity marker (through *UAS-GFP* expression, for example) or is a tool to drive other *UAS-transgenes* (eg., *dsRNA*) for genetic analysis; simultaneous use can make phenotypic interpretations difficult. Another key issue common to the use of *Gal4*-based drivers in *Drosophila* is that it is nearly impossible to perform and/or interpret results from experiments that combine multiple *Gal4* reporters into the same background, as expression patterns will combine. In short, continued progress in characterizing the fly hematopoietic system will require the implementation of new or modified genetic technologies that can complement existing tools. In the past few years, several groups have expanded the *Drosophila* molecular genetic tool kit through the generation of new bipartite gene expression systems [103], such as the *LexA-lexAop* [104] and *QF-QUAS* systems [105]. Once adapted for use in the hematopoietic system (through the creation of appropriate driver lines such as *Hml-QF* or *dome-QF*), it will be possible to use such lines in combination with the *Gal4/UAS* system to enhance and refine analyses.

An alternative approach taken by our lab to overcome current limitations and to facilitate genetic analysis in the lymph gland has been to generate a new line of flies that uses multiple *Gal4*-independent fluorescent protein reporters to simultaneously monitor different lymph gland cell populations. To generate this new hematopoietic tool, we first re-tasked the

*dome-MESO* progenitor cell enhancer [83] to drive expression of enhanced blue fluorescent protein 2 (EBFP2, Ex/Em: 383/448 nm) [106]. This transgene was then recombined with differentiating-cell marker *Hemolectin-DsRed* [58] and the PSC/niche cell marker *hedgehog-GFP* [81] onto the X chromosome. As can be seen in Figure 1, each reporter clearly identifies their respective cell in the lymph gland and, therefore, will be extremely useful for assessing phenotypes in various experimental contexts, such as genetic screening. This triple recombinant marker line, as well as the new, individual fluorescent *dome-MESO* transgenic lines, will be made freely available to the research community.

In contrast to other organs such as the imaginal discs, brain, gut, and fat body, the lymph gland is likely to be unfamiliar to most *Drosophilists*. Because of its relatively small size and delicate nature, the lymph gland presents a significant challenge to anyone interested in analyzing it by dissection and microscopy. To provide guidance in this area, a protocol for lymph gland dissection is described, and is followed by a standard immunostaining procedure and a summary of how to mount high quality, intact lymph glands on glass slides for microscopic analysis (see Boxes 2, 3, and 4).

## Analysis of adult hemocytes

As mentioned, relatively little is known about adult hemocytes or hematopoiesis, owing primarily to a lack of useful tools and methods. In contrast to larvae, which are akin to fluid-filled balloons, adult bodies are rigid and have a relatively lower hemolymph volume, making analysis of circulating cells by “bleeding” more challenging. This becomes even more difficult in light of the fact that many, if not most, adult hemocytes are sessile, associated with the adult body wall and internal organs and structures. Many such adult hemocytes can be visualized through the cuticle in whole-mount preparations using markers such as *Hml-gal4 UAS-GFP*, however approach this is not very amenable to more detailed analyses of blood cell fate, proliferation, or function. It is clear that more work needs to be done to develop reagents and techniques to make the adult blood system more accessible to molecular genetic analysis. To date, the most commonly applied method for visualizing and analyzing adult hemocytes is to perfuse the adult hemocoel with aqueous buffer, which flushes hemocytes out of the adult and onto a slide for analysis. An adult perfusion protocol is described in Box 7. Recently, such a method was used in the identification of the first adult-specific blood marker, Ad1 [107], and further demonstrated that adult hemocytes maintain NimC1 (P1) expression but lose expression of Hemese (He), the blood-specific pan hemocyte marker found in larval hemocytes.

## Summary

As with any developing field of study, continued effort on the part of researchers will enhance the ability to manipulate and analyze the fly hematopoietic system. This advancement will undoubtedly rely in large part on the identification of improved molecular markers offering higher resolution identification of cell types within the hematopoietic lineages, as well as on improved technologies for genetic analysis and imaging. Of particular importance will be approaches that allow for analyses *in vivo* during the larval, pupal, and adult stages, which may mitigate current difficulties in extracting blood cells from animals

for *ex vivo* analyses. This review has highlighted some of the more useful genetic markers and analytic methods, both by pointing toward expertise in the published literature where appropriate and by providing here a set of basic protocols for dealing with the larval and adult blood system.

## Acknowledgments

Because of space limitations, the review could not be exhaustive in its description of the hematopoietic system. We apologize to those whose important work in the field is not cited here. This work was supported by NIH grants R01 HL067395 and R21 AI094048 to U.B., and a China Scholarship Council award to T.L.

## References

1. Rizki, TM. Blood cells of *Drosophila* as related to metamorphosis.. In: Campbell, FL., editor. *Physiology of Insect Development*. Chicago University Press; Chicago, IL.: 1956. p. 91-94.
2. Rizki TM, Rizki RM. Properties of the larval hemocytes of *Drosophila melanogaster*. *Experientia*. 1980; 36:1223–1226.
3. Till JE, McCulloch EA. A direct measurement of the radiation sensitivity of normal mouse bone marrow cells. *Radiation Research*. 1961; 14:213–222. [PubMed: 13776896]
4. Rizki, TM. The circulatory system and associated cells and tissues. In: Ashburner, M.; Wright, TRF., editors. *The Genetics and Biology of Drosophila*. Academic Press; New York, London: 1978.
5. Lanot R, Zachary D, Holder F, Meister M. Postembryonic hematopoiesis in *Drosophila*. *Dev Biol*. 2001; 230(2):243–57. [PubMed: 11161576]
6. Rizki TM, Rizki RM. Lamellocyte differentiation in *Drosophila* larvae parasitized by *Leptopilina*. *Dev Comp Immunol*. 1992; 16(2-3):103–10. [PubMed: 1499832]
7. Sorrentino RP, Carton Y, Govind S. Cellular immune response to parasite infection in the *Drosophila* lymph gland is developmentally regulated. *Dev Biol*. 2002; 243(1):65–80. [PubMed: 11846478]
8. Evans CJ, Hartenstein V, Banerjee U. Thicker than blood: conserved mechanisms in *Drosophila* and vertebrate hematopoiesis. *Dev Cell*. 2003; 5(5):673–90. [PubMed: 14602069]
9. Gupta, AP. Hemocyte types: their structures, synonymies, interrelationships, and taxonomic significance.. In: Gupta, AP., editor. *Insect Hemocytes: Development, forms, functions, and techniques*. Cambridge University Press; 1979. p. 85-127.
10. Shrestha R, Gateff E. Ultrastructure and cytochemistry of the cell types of the larval hematopoietic organs and hemolymph of *Drosophila melanogaster*. *Dev. Growth Differ*. 1982; 24:65–82.
11. Agaisse H, Petersen UM, Boutros M, Mathey-Prevot B, Perrimon N. Signaling role of hemocytes in *Drosophila* JAK/STAT-dependent response to septic injury. *Dev Cell*. 2003; 5(3):441–50. [PubMed: 12967563]
12. Hetru C, Troxler L, Hoffmann JA. *Drosophila melanogaster* antimicrobial defense. *J Infect Dis*. 2003; 187(Suppl 2):S327–34. [PubMed: 12792847]
13. Nappi AJ. Inhibition by parasites of melanotic tumour formation in *Drosophila melanogaster*. *Nature*. 1975; 255(5507):402–4. [PubMed: 805380]
14. Samakovlis C, Kimbrell DA, Kylsten P, Engstrom A, Hultmark D. The immune response in *Drosophila*: pattern of cecropin expression and biological activity. *EMBO J*. 1990; 9(9):2969–76. [PubMed: 2390977]
15. Tzou P, Reichhart JM, Lemaitre B. Constitutive expression of a single antimicrobial peptide can restore wild-type resistance to infection in immunodeficient *Drosophila* mutants. *Proc Natl Acad Sci U S A*. 2002; 99(4):2152–7. [PubMed: 11854512]
16. Brennan CA, Delaney JR, Schneider DS, Anderson KV. Psidin is required in *Drosophila* blood cells for both phagocytic degradation and immune activation of the fat body. *Curr Biol*. 2007; 17(1):67–72. [PubMed: 17208189]
17. Elrod-Erickson M, Mishra S, Schneider D. Interactions between the cellular and humoral immune responses in *Drosophila*. *Curr Biol*. 2000; 10(13):781–4. [PubMed: 10898983]

18. Matova N, Anderson KV. Rel/NF-kappaB double mutants reveal that cellular immunity is central to *Drosophila* host defense. *Proc Natl Acad Sci U S A*. 2006; 103(44):16424–9. [PubMed: 17060622]
19. Pham LN, Dionne MS, Shirasu-Hiza M, Schneider DS. A specific primed immune response in *Drosophila* is dependent on phagocytes. *PLoS Pathog*. 2007; 3(3):e26. [PubMed: 17352533]
20. Fessler LI, Nelson RE, Fessler JH. *Drosophila* extracellular matrix. *Methods Enzymol*. 1994; 245:271–94. [PubMed: 7760738]
21. Franc NC. Phagocytosis of apoptotic cells in mammals, *Caenorhabditis elegans* and *Drosophila melanogaster*: molecular mechanisms and physiological consequences. *Front Biosci*. 2002; 7:d1298–313. [PubMed: 11991836]
22. Murray MA, Fessler LI, Palka J. Changing distributions of extracellular matrix components during early wing morphogenesis in *Drosophila*. *Dev Biol*. 1995; 168(1):150–65. [PubMed: 7883070]
23. Yasothornsrikul S, Davis WJ, Cramer G, Kimbrell DA, Dearolf CR. viking: identification and characterization of a second type IV collagen in *Drosophila*. *Gene*. 1997; 198(1-2):17–25. [PubMed: 9370260]
24. Sears HC, Kennedy CJ, Garrity PA. Macrophage-mediated corpse engulfment is required for normal *Drosophila* CNS morphogenesis. *Development*. 2003; 130(15):3557–65. [PubMed: 12810602]
25. Franc NC, Dimarcq JL, Lagueux M, Hoffmann J, Ezekowitz RA. Croquemort, a novel *Drosophila* hemocyte/macrophage receptor that recognizes apoptotic cells. *Immunity*. 1996; 4(5):431–43. [PubMed: 8630729]
26. Fadok VA, Bratton DL, Rose DM, Pearson A, Ezekowitz RA, et al. A receptor for phosphatidylserine-specific clearance of apoptotic cells. *Nature*. 2000; 405(6782):85–90. [PubMed: 10811223]
27. Freeman MR, Delrow J, Kim J, Johnson E, Doe CQ. Unwrapping glial biology: Gcm target genes regulating glial development, diversification, and function. *Neuron*. 2003; 38(4):567–80. [PubMed: 12765609]
28. Manaka J, Kuraishi T, Shiratsuchi A, Nakai Y, Higashida H, et al. Draper-mediated and phosphatidylserine-independent phagocytosis of apoptotic cells by *Drosophila* hemocytes/macrophages. *J Biol Chem*. 2004; 279(46):48466–76. [PubMed: 15342648]
29. Ramet M, Manfrulli P, Pearson A, Mathey-Prevot B, Ezekowitz RA. Functional genomic analysis of phagocytosis and identification of a *Drosophila* receptor for *E. coli*. *Nature*. 2002; 416(6881):644–8. [PubMed: 11912489]
30. Ramet M, Pearson A, Manfrulli P, Li X, Koziel H, et al. *Drosophila* scavenger receptor CI is a pattern recognition receptor for bacteria. *Immunity*. 2001; 15(6):1027–38. [PubMed: 11754822]
31. Royet J, Reichhart JM, Hoffmann JA. Sensing and signaling during infection in *Drosophila*. *Curr Opin Immunol*. 2005; 17(1):11–7. [PubMed: 15653304]
32. Kocks C, Cho JH, Nehme N, Ulvila J, Pearson AM, et al. Eater, a transmembrane protein mediating phagocytosis of bacterial pathogens in *Drosophila*. *Cell*. 2005; 123(2):335–46. [PubMed: 16239149]
33. Kurucz E, Markus R, Zsomboki J, Folkl-Medzihradsky K, Darula Z, et al. Nimrod, a putative phagocytosis receptor with EGF repeats in *Drosophila* plasmatocytes. *Curr Biol*. 2007; 17(7):649–54. [PubMed: 17363253]
34. Watson FL, Puttmann-Holgado R, Thomas F, Lamar DL, Hughes M, et al. Extensive diversity of Ig-superfamily proteins in the immune system of insects. *Science*. 2005; 309(5742):1874–8. [PubMed: 16109846]
35. Schmucker D, Clemens JC, Shu H, Worby CA, Xiao J, et al. *Drosophila* Dscam is an axon guidance receptor exhibiting extraordinary molecular diversity. *Cell*. 2000; 101(6):671–84. [PubMed: 10892653]
36. Castillejo-Lopez C, Hacker U. The serine protease Sp7 is expressed in blood cells and regulates the melanization reaction in *Drosophila*. *Biochem Biophys Res Commun*. 2005; 338(2):1075–82. [PubMed: 16256951]

37. Chosa N, Fukumitsu T, Fujimoto K, Ohnishi E. Activation of prophenoloxidase A1 by an activating enzyme in *Drosophila melanogaster*. *Insect Biochem Mol Biol*. 1997; 27(1):61–8. [PubMed: 9061929]
38. De Gregorio E, Spellman PT, Tzou P, Rubin GM, Lemaitre B. The Toll and Imd pathways are the major regulators of the immune response in *Drosophila*. *EMBO J*. 2002; 21(11):2568–79. [PubMed: 12032070]
39. Tang H, Kambris Z, Lemaitre B, Hashimoto C. Two proteases defining a melanization cascade in the immune system of *Drosophila*. *J Biol Chem*. 2006; 281(38):28097–104. [PubMed: 16861233]
40. Rizki TM, Rizki RM. Paracrystalline inclusions of *D. melanogaster* hemocytes have prophenoloxidases. *Genetics*. 1985; 110(S98)
41. Bidla G, Dushay MS, Theopold U. Crystal cell rupture after injury in *Drosophila* requires the JNK pathway, small GTPases and the TNF homolog Eiger. *J Cell Sci*. 2007; 120(Pt 7):1209–15. [PubMed: 17356067]
42. Peebles EE, Geisler A, Whitcraft CJ, Oliver CP. Comparative studies of phenol oxidase activity during pupal development of three lozenge mutants (lz8,lz,lzk) of *Drosophila melanogaster*. *Genetics*. 1969; 62(1):161–70. [PubMed: 4984328]
43. Rizki TM, Rizki RM. Alleles of lz as suppressors of the Bc-phenix in *Drosophila melanogaster*. *Genetics*. 1981; 97(S90)
44. Bidla G, Lindgren M, Theopold U, Dushay MS. Hemolymph coagulation and phenoloxidase in *Drosophila* larvae. *Dev Comp Immunol*. 2005; 29(8):669–79. [PubMed: 15854679]
45. Russo J, Dupas S, Frey F, Carton Y, Brehelin M. Insect immunity: early events in the encapsulation process of parasitoid (*Leptopilina boulardi*) eggs in resistant and susceptible strains of *Drosophila*. *Parasitology*. 1996; 112(Pt 1):135–42. [PubMed: 8587797]
46. Luo H, Rose PE, Roberts TM, Dearolf CR. The Hopscotch Jak kinase requires the Raf pathway to promote blood cell activation and differentiation in *Drosophila*. *Mol Genet Genomics*. 2002; 267(1):57–63. [PubMed: 11919715]
47. Markus R, Kurucz E, Rus F, Ando I. Sterile wounding is a minimal and sufficient trigger for a cellular immune response in *Drosophila melanogaster*. *Immunol Lett*. 2005; 101(1):108–11. [PubMed: 15964636]
48. Lemaitre B. The road to Toll. *Nat Rev Immunol*. 2004; 4(7):521–7. [PubMed: 15229471]
49. Romeo Y, Lemaitre B. *Drosophila* immunity: methods for monitoring the activity of Toll and Imd signaling pathways. *Methods Mol Biol*. 2008; 415:379–94. [PubMed: 18370166]
50. Small C, Paddibhatla I, Rajwani R, Govind S. An introduction to parasitic wasps of *Drosophila* and the antiparasite immune response. *J Vis Exp*. 2012; (63):e3347. [PubMed: 22588641]
51. Hartenstein V. Blood cells and blood cell development in the animal kingdom. *Annu Rev Cell Dev Biol*. 2006; 22:677–712. [PubMed: 16824014]
52. Lebestky T, Chang T, Hartenstein V, Banerjee U. Specification of *Drosophila* hematopoietic lineage by conserved transcription factors. *Science*. 2000; 288(5463):146–9. [PubMed: 10753120]
53. Tepass U, Fessler LI, Aziz A, Hartenstein V. Embryonic origin of hemocytes and their relationship to cell death in *Drosophila*. *Development*. 1994; 120(7):1829–37. [PubMed: 7924990]
54. Jung SH, Evans CJ, Uemura C, Banerjee U. The *Drosophila* lymph gland as a developmental model of hematopoiesis. *Development*. 2005; 132(11):2521–33. [PubMed: 15857916]
55. Rugendorff AE, Younossi-Hartenstein A, Hartenstein V. Embryonic origin and differentiation of the *Drosophila* heart. *Roux Arch. Dev. Biol*. 1994; 203:266–280.
56. Holz A, Bossinger B, Strasser T, Janning W, Klapper R. The two origins of hemocytes in *Drosophila*. *Development*. 2003; 130(20):4955–62. [PubMed: 12930778]
57. Robertson CW. The metamorphosis of *Drosophila melanogaster*, including an accurately timed account of the principal morphological changes. *J. Morphol*. 1936; 59:351–399.
58. Makhijani K, Alexander B, Tanaka T, Rulifson E, Bruckner K. The peripheral nervous system supports blood cell homing and survival in the *Drosophila* larva. *Development*. 2011; 138(24):5379–91. [PubMed: 22071105]

59. Markus R, Laurinyecz B, Kurucz E, Honti V, Bajusz I, et al. Sessile hemocytes as a hematopoietic compartment in *Drosophila melanogaster*. *Proc Natl Acad Sci U S A*. 2009; 106(12):4805–9. [PubMed: 19261847]
60. Stofanko M, Kwon SY, Badenhorst P. A misexpression screen to identify regulators of *Drosophila* larval hemocyte development. *Genetics*. 2008; 180(1):253–67. [PubMed: 18757933]
61. Waltzer L, Ferjoux G, Bataille L, Haenlin M. Cooperation between the GATA and RUNX factors Serpent and Lozenge during *Drosophila* hematopoiesis. *EMBO J*. 2003; 22(24):6516–25. [PubMed: 14657024]
62. Milchanowski AB, Henkenius AL, Narayanan M, Hartenstein V, Banerjee U. Identification and characterization of genes involved in embryonic crystal cell formation during *Drosophila* hematopoiesis. *Genetics*. 2004; 168(1):325–39. [PubMed: 15454546]
63. Crozatier M, Ubeda JM, Vincent A, Meister M. Cellular immune response to parasitization in *Drosophila* requires the EBF orthologue collier. *PLoS Biol*. 2004; 2(8):E196. [PubMed: 15314643]
64. Mandal L, Martinez-Agosto JA, Evans CJ, Hartenstein V, Banerjee U. A Hedgehog- and Antennapedia-dependent niche maintains *Drosophila* haematopoietic precursors. *Nature*. 2007; 446(7133):320–4. [PubMed: 17361183]
65. Evans IR, Hu N, Skaer H, Wood W. Interdependence of macrophage migration and ventral nerve cord development in *Drosophila* embryos. *Development*. 2010; 137(10):1625–33. [PubMed: 20392742]
66. Tucker PK, Evans IR, Wood W. Ena drives invasive macrophage migration in *Drosophila* embryos. *Dis Model Mech*. 2011; 4(1):126–34. [PubMed: 21045209]
67. Wood W, Faria C, Jacinto A. Distinct mechanisms regulate hemocyte chemotaxis during development and wound healing in *Drosophila melanogaster*. *J Cell Biol*. 2006; 173(3):405–16. [PubMed: 16651377]
68. Stramer B, Wood W, Galko MJ, Redd MJ, Jacinto A, et al. Live imaging of wound inflammation in *Drosophila* embryos reveals key roles for small GTPases during in vivo cell migration. *J Cell Biol*. 2005; 168(4):567–73. [PubMed: 15699212]
69. Evans IR, Zanet J, Wood W, Stramer BM. Live imaging of *Drosophila melanogaster* embryonic hemocyte migrations. *J Vis Exp*. 2010; (36)
70. Moreira CG, Regan JC, Zaidman-Remy A, Jacinto A, Prag S. *Drosophila* hemocyte migration: an in vivo assay for directional cell migration. *Methods Mol Biol*. 2011; 769:249–60. [PubMed: 21748681]
71. Sampson CJ, Williams MJ. Protocol for ex vivo incubation of *Drosophila* primary post-embryonic haemocytes for real-time analyses. *Methods Mol Biol*. 2012; 827:359–67. [PubMed: 22144286]
72. Sinenko SA, Hung T, Moroz T, Tran QM, Sidhu S, et al. Genetic manipulation of AML1-ETO-induced expansion of hematopoietic precursors in a *Drosophila* model. *Blood*. 2010; 116(22):4612–20. [PubMed: 20688956]
73. Moreira CG, Jacinto A, Prag S. *Drosophila* integrin adhesion complexes are essential for hemocyte migration in vivo. *Biol Open*. 2013; 2(8):795–801. [PubMed: 23951405]
74. Tirouvanziam R, Davidson CJ, Lipsick JS, Herzenberg LA. Fluorescence-activated cell sorting (FACS) of *Drosophila* hemocytes reveals important functional similarities to mammalian leukocytes. *Proc Natl Acad Sci U S A*. 2004; 101(9):2912–7. [PubMed: 14976247]
75. Vilmos P, Nagy I, Kurucz E, Hultmark D, Gateff E, et al. A rapid rosetting method for separation of hemocyte sub-populations of *Drosophila melanogaster*. *Dev Comp Immunol*. 2004; 28(6):555–63. [PubMed: 15177109]
76. Zettervall CJ, Anderl I, Williams MJ, Palmer R, Kurucz E, et al. A directed screen for genes involved in *Drosophila* blood cell activation. *Proc Natl Acad Sci U S A*. 2004; 101(39):14192–7. [PubMed: 15381778]
77. Crozatier M, Vincent A. *Drosophila*: a model for studying genetic and molecular aspects of haematopoiesis and associated leukaemias. *Dis Model Mech*. 2011; 4(4):439–45. [PubMed: 21669932]

78. Krzemien J, Dubois L, Makki R, Meister M, Vincent A, et al. Control of blood cell homeostasis in *Drosophila* larvae by the posterior signalling centre. *Nature*. 2007; 446(7133):325–8. [PubMed: 17361184]
79. Lebestky T, Jung SH, Banerjee U. A Serrate-expressing signaling center controls *Drosophila* hematopoiesis. *Genes Dev*. 2003; 17(3):348–53. [PubMed: 12569125]
80. Pennetier D, Oyallon J, Morin-Poulard I, Dejean S, Vincent A, et al. Size control of the *Drosophila* hematopoietic niche by bone morphogenetic protein signaling reveals parallels with mammals. *Proc Natl Acad Sci U S A*. 2012; 109(9):3389–94. [PubMed: 22331866]
81. Tokusumi T, Shoue DA, Tokusumi Y, Stoller JR, Schulz RA. New hemocyte-specific enhancer-reporter transgenes for the analysis of hematopoiesis in *Drosophila*. *Genesis*. 2009; 47(11):771–4. [PubMed: 19830816]
82. Bourbon HM, Gonzy-Treboul G, Peronnet F, Alin MF, Ardourel C, et al. A P-insertion screen identifying novel X-linked essential genes in *Drosophila*. *Mech Dev*. 2002; 110(1-2):71–83. [PubMed: 11744370]
83. Hombria JC, Brown S, Hader S, Zeidler MP. Characterisation of Upd2, a *Drosophila* JAK/STAT pathway ligand. *Dev Biol*. 2005; 288(2):420–33. [PubMed: 16277982]
84. Mondal BC, Mukherjee T, Mandal L, Evans CJ, Sinenko SA, et al. Interaction between differentiating cell- and niche-derived signals in hematopoietic progenitor maintenance. *Cell*. 2011; 147(7):1589–600. [PubMed: 22196733]
85. Sinenko SA, Mandal L, Martinez-Agosto JA, Banerjee U. Dual role of wingless signaling in stem-like hematopoietic precursor maintenance in *Drosophila*. *Dev Cell*. 2009; 16(5):756–63. [PubMed: 19460351]
86. Avet-Rochex A, Boyer K, Polesello C, Gobert V, Osman D, et al. An in vivo RNA interference screen identifies gene networks controlling *Drosophila melanogaster* blood cell homeostasis. *BMC Dev Biol*. 2010; 10:65. [PubMed: 20540764]
87. Shim J, Mukherjee T, Banerjee U. Direct sensing of systemic and nutritional signals by haematopoietic progenitors in *Drosophila*. *Nat Cell Biol*. 2012; 14(4):394–400. [PubMed: 22407365]
88. Tokusumi T, Tokusumi Y, Hopkins DW, Shoue DA, Corona L, et al. Germ line differentiation factor Bag of Marbles is a regulator of hematopoietic progenitor maintenance during *Drosophila* hematopoiesis. *Development*. 2011; 138(18):3879–84. [PubMed: 21813570]
89. Owusu-Ansah E, Banerjee U. Reactive oxygen species prime *Drosophila* haematopoietic progenitors for differentiation. *Nature*. 2009; 461(7263):537–41. [PubMed: 19727075]
90. Robinson KM, Janes MS, Pehar M, Monette JS, Ross MF, et al. Selective fluorescent imaging of superoxide in vivo using ethidium-based probes. *Proc Natl Acad Sci U S A*. 2006; 103(41):15038–43. [PubMed: 17015830]
91. Owusu-Ansah, E.; Yavari, A.; Banerjee, U. A protocol for \_in vivo\_ detection of reactive oxygen species. 2008.
92. Sinenko SA, Mathey-Prevot B. Increased expression of *Drosophila* tetraspanin, Tsp68C, suppresses the abnormal proliferation of ytr-deficient and Ras/Raf-activated hemocytes. *Oncogene*. 2004; 23(56):9120–8. [PubMed: 15480416]
93. Galko MJ, Krasnow MA. Cellular and genetic analysis of wound healing in *Drosophila* larvae. *PLoS Biol*. 2004; 2(8):E239. [PubMed: 15269788]
94. Nelson RE, Fessler LI, Takagi Y, Blumberg B, Keene DR, et al. Peroxidase: a novel enzyme-matrix protein of *Drosophila* development. *EMBO J*. 1994; 13(15):3438–47. [PubMed: 8062820]
95. Asha H, Nagy I, Kovacs G, Stetson D, Ando I, et al. Analysis of Ras-induced overproliferation in *Drosophila* hemocytes. *Genetics*. 2003; 163(1):203–15. [PubMed: 12586708]
96. Kurucz E, Vaczi B, Markus R, Laurinyecz B, Vilmos P, et al. Definition of *Drosophila* hemocyte subsets by cell-type specific antigens. *Acta Biol Hung*. 2007; 58(Suppl):95–111. [PubMed: 18297797]
97. Honti V, Cinege G, Csordas G, Kurucz E, Zsamboki J, et al. Variation of NimC1 expression in *Drosophila* stocks and transgenic strains. *Fly (Austin)*. 2013; 7(4)

98. Braun A, Lemaitre B, Lanot R, Zachary D, Meister M. *Drosophila* immunity: analysis of larval hemocytes by P-element-mediated enhancer trap. *Genetics*. 1997; 147(2):623–34. [PubMed: 9335599]
99. Tokusumi T, Sorrentino RP, Russell M, Ferrarese R, Govind S, et al. Characterization of a lamellocyte transcriptional enhancer located within the *misshapen* gene of *Drosophila melanogaster*. *PLoS One*. 2009; 4(7):e6429. [PubMed: 19641625]
100. Rus F, Kurucz E, Markus R, Sinenko SA, Laurinyecz B, et al. Expression pattern of Filamin-240 in *Drosophila* blood cells. *Gene Expr Patterns*. 2006; 6(8):928–34. [PubMed: 16616709]
101. Irving P, Ubeda JM, Doucet D, Troxler L, Lagueux M, et al. New insights into *Drosophila* larval haemocyte functions through genome-wide analysis. *Cell Microbiol*. 2005; 7(3):335–50. [PubMed: 15679837]
102. Honti V, Kurucz E, Csordas G, Laurinyecz B, Markus R, et al. In vivo detection of lamellocytes in *Drosophila melanogaster*. *Immunol Lett*. 2009; 126(1-2):83–4. [PubMed: 19695290]
103. del Valle Rodriguez A, Didiano D, Desplan C. Power tools for gene expression and clonal analysis in *Drosophila*. *Nat Methods*. 2012; 9(1):47–55. [PubMed: 22205518]
104. Lai SL, Lee T. Genetic mosaic with dual binary transcriptional systems in *Drosophila*. *Nat Neurosci*. 2006; 9(5):703–9. [PubMed: 16582903]
105. Potter CJ, Tasic B, Russler EV, Liang L, Luo L. The Q system: a repressible binary system for transgene expression, lineage tracing, and mosaic analysis. *Cell*. 2010; 141(3):536–48. [PubMed: 20434990]
106. Ai HW, Shaner NC, Cheng Z, Tsien RY, Campbell RE. Exploration of new chromophore structures leads to the identification of improved blue fluorescent proteins. *Biochemistry*. 2007; 46(20):5904–10. [PubMed: 17444659]
107. Honti V, Csordas G, Kurucz E, Markus R, Ando I. The cell-mediated immunity of *Drosophila melanogaster*: hemocyte lineages, immune compartments, microanatomy and regulation. *Dev Comp Immunol*. 2014; 42(1):47–56. [PubMed: 23800719]
108. Bruckner K, Kockel L, Duchek P, Luque CM, Rorth P, et al. The PDGF/VEGF receptor controls blood cell survival in *Drosophila*. *Dev Cell*. 2004; 7(1):73–84. [PubMed: 15239955]
109. Bernardoni R, Vivancos V, Giangrande A. *glide/gcm* is expressed and required in the scavenger cell lineage. *Dev. Biol*. 1997; 191:118–130. [PubMed: 9356176]
110. Cho NK, Keyes L, Johnson E, Heller J, Ryner L, et al. Developmental control of blood cell migration by the *Drosophila* VEGF pathway. 108. *Cell*. 2002:865–876.
111. Olofsson B, Page DT. Condensation of the central nervous system in embryonic *Drosophila* is inhibited by blocking hemocyte migration or neural activity. *Dev Biol*. 2005; 279(1):233–243. [PubMed: 15708571]
112. Waltzer L, Bataillé L, Peyrefitte S, Haenlin M. Two isoforms of *Serpent* containing either one or two GATA zinc fingers have different roles in *Drosophila* haematopoiesis. *EMBO J*. 2002; 21(20):5477–86. [PubMed: 12374748]
113. Ward EJ, Skeath JB. Characterization of a novel subset of cardiac cells and their progenitors in the *Drosophila* embryo. *Development*. 2000; 127(22):4959–69. [PubMed: 11044409]
114. Han Z, Olson EN. Hand is a direct target of Tinman and GATA factors during *Drosophila* cardiogenesis and hematopoiesis. *Development*. 2005; 132(15):3525–36. [PubMed: 15975941]
115. Evans CJ, Olson JM, Ngo KT, Kim E, Lee NE, et al. G-TRACE: rapid Gal4-based cell lineage analysis in *Drosophila*. *Nat Methods*. 2009; 6(8):603–5. [PubMed: 19633663]
116. Kimbrell DA, Hice C, Bolduc C, Kleinhesselink K, Beckingham K. The Dorothy enhancer has Tinman binding sites and drives hopscotch-induced tumor formation. *Genesis*. 2002; 34:23–28. [PubMed: 12324942]
117. Kurucz E, Zettervall CJ, Sinka R, Vilmos P, Pivarcsi A, et al. Hemese, a hemocyte-specific transmembrane protein, affects the cellular immune response in *Drosophila*. *Proc Natl Acad Sci U S A*. 2003; 100(5):2622–7. [PubMed: 12598653]
118. Shim J, Mukherjee T, Mondal BC, Liu T, Young GC, et al. Olfactory control of blood progenitor maintenance. *Cell*. 2013; 155(5):1141–53. [PubMed: 24267893]

119. Sinenko SA, Kim EK, Wynn R, Manfrulli P, Ando I, et al. Yantar, a conserved arginine-rich protein is involved in *Drosophila* hemocyte development. *Dev Biol.* 2004; 273(1):48–62. [PubMed: 15302597]
120. Rizki TM, Rizki RM, Grell EH. A mutant affecting the crystal cells of *Drosophila melanogaster*. *Roux Arch. Dev. Biol.* 1980; 188:91–99.
121. Mukherjee T, Kim WS, Mandal L, Banerjee U. Interaction between Notch and Hif-alpha in development and survival of *Drosophila* blood cells. *Science.* 2011; 332(6034):1210–3. [PubMed: 21636775]

**Box 1 | Collection and processing of larval circulating blood cells**

1. Collect larvae and wash thoroughly. This is important to ensure that the circulating cell preparation is not contaminated with food debris, yeast, and bacteria (which will obscure observation of the cells) from the culture.
2. Place a 14-well glass slide on the stereo microscope and illuminate with the transmitted light base.
3. Place a 20  $\mu$ L drop of 1XPBS on one well spot of a 14-well glass slide. Do not place a drop on each well at the beginning as evaporation will significantly decrease the buffer volume as you proceed. Instead, use one drop at a time.
4. Place one larva in the drop of buffer.
5. Using two pairs of forceps, pinch the body wall/cuticle at the posterior of the larva. Keep one forceps in place while tearing the cuticle down the length of the larva toward the anterior. As the hemocoel is opened, hemolymph will begin to stream into the 1XPBS buffer; this should be relatively easy to observe because of differential fluid densities. To obtain a clean preparation of isolated blood cells, it is essential to grab and tear only the cuticle. Avoid damaging internal organs particularly the gut, the contents of which will significantly contaminate the sample and will be difficult to wash away from the slide surface.
6. Once the larval hemocoel is open, use forceps to gently move the larva through the buffer to evenly distribute the hemolymph and hemocytes throughout.
7. Using single tines from two different forceps, gently remove the carcass from the drop of buffer containing the blood cells. Do this slowly so that minimal buffer is removed. Do NOT use one pair of forceps to pick up the carcass as too much buffer will also be removed due to capillary action between the tines.
8. Repeat for remaining samples.
9. Carefully place slide into humidified chamber. Incubate for 30 minutes to let blood cells settle onto the glass surface.
10. Remove slide from the humidified chamber.
11. Under the microscope, carefully remove the 1XPBS buffer from each drop using a pipette.
12. Add 20  $\mu$ L of freshly prepared fixative to each well and place into humidified chamber for 10 minutes at room temperature.
13. Under the microscope, remove fixative and replace with 20  $\mu$ L of 1 $\times$ PBS for short term storage at 4C or proceed with subsequent assay of choice, such as immunostaining. For steps using buffers with detergents, volumes of 5-10  $\mu$ L are recommended to avoid mixing between wells due to loss of surface tension. All steps should be carried out in a humidified chamber to avoid evaporation.

14. After processing, place a small drop of VECTASHIELD on each well, followed by a glass cover slip.
15. Seal with nail polish around the edges and proceed to imaging.

### Box 2 | Lymph gland dissection

1. Using a large transfer pipette, place several large drops of 1XPBS on the silicone dissection plate.
2. Place one larva in each drop. As technique improves, several larvae can be placed in each drop without interference.
3. Place the plate on a standard stereo dissecting microscope with a transmitted light base.
4. Illuminate the sample. If lighting makes use of standard light bulbs (instead of a cold LED source), take care to work quickly so that the larvae/samples do not overheat. Dissect for a maximum of 30 minutes, with samples stored on ice, before fixation.
5. Using forceps, orient larva dorsal-side up, anterior to the right (for a right-handed person).
6. Using the left hand, gently hold larvae at 75% length (A=0%, P=100%).
7. Using the right hand, pinch the dorsal cuticle (only, no internal organs) at 70-75% length, and slowly pull toward the anterior. Pull the entire length, stopping at the mouth hooks. This will open the body cavity without damaging internal structures. Sometimes the dorsal cuticle will completely detach or the ventral cuticle also; either is ok as long as the body cavity is open without disrupting internal organs. *Alternative method:* Using both pairs of forceps, rip the larva in half, then invert the anterior part by pushing the mouth hooks into the body cavity longitudinally (imagine inverting a sock by pushing your hand inward at the toes). Once inverted, proceed with step 10 below. Note: this method opens the gut, which can make subsequent steps more difficult because food debris, yeast, and bacteria cloud the buffer.
8. Use both forceps to sever the larvae (internal organs and ventral cuticle) completely at 75% length (near where the larva was being held by the left hand).
9. Stabilize the anterior portion of the larva with the left hand by grabbing the ventral cuticle. Using the right hand forceps, reach inside the body cavity and directly grasp the mouth hooks (do not grab them externally as this will include cuticle).
10. Gently pull the mouth hooks out, separating them from the body wall and cuticle. Removing the mouth hooks in this manner will bring along a complex of structures that includes the eye imaginal discs, the brain, the dorsal vessel, the prothoracic (ring) gland, and the lymph gland; salivary glands are also commonly included.
11. Gently separate these lymph gland-containing complexes from any residual tissues such as fat body.

- 12.** Carefully remove the complex in a small drop of 1XPBS retained between the tines of a single forceps) to a microfuge tube containing cold 1mL 1XPBS + 1 drop 1XPBST on ice. The PBST break the surface tension of the buffer, allowing samples to sink to the bottom of the tube.
- 13.** Collect as many lymph gland complexes as needed.
- 14.** Replace buffer with 1 mL of freshly prepared fixative.
- 15.** Place microfuge tube on mixer (Nutator, orbital shaker, or similar) for 30 minutes.
- 16.** Proceed to immunostaining (described below) or other analytical method.
- 17.** For short term storage, wash fixed tissue with 1 mL 1XPBS for 5 minutes; repeat; store at 4°C.

### Box 3 | Lymph gland immunostaining

For this procedure, lymph glands are left attached to the mouth hook/eye disc/brain complex, 1) to provide a “handle” so that lymph glands are not damaged because of physical manipulation and 2) because these tissues are often convenient immunostaining controls

1. With lymph glands either in fix or 1XPBS in a 1.5 mL microfuge tube, discard buffer.
2. Add 1mL 1XPBST (1XPBS + 0.4% Triton), place on mixer (Nutator or similar) for 15 minutes; repeat twice.
3. During the washing stage, prepare 10% normal goat serum (NGS)/1xPBST blocking solution.
4. Discard wash buffer from lymph glands and replace with 1 mL blocking solution; block in for 30 minutes.
5. Dilute primary antibody in blocking solution at desired concentration.
6. Discard blocking solution from lymph glands and replace with 100  $\mu$ L or greater volume of primary antibody and mix well. For smaller volumes and/or limited primary antibody, 0.6 mL microfuge tubes, PCR tubes, or Terasaki microtiter plates (with 10  $\mu$ L wells) can be used.
7. Place tubes in a rack. For short term incubations at room temperature, intermittent mixing or placement on an orbital shaker is helpful.
8. Incubate for 3 hours at room temperature or overnight at 4C.
9. Remove primary antibody.
10. Add 1mL 1XPBST, place on mixer for 15 minutes; repeat twice.
11. Discard wash buffer from lymph glands and replace with 1 mL blocking solution; block 10-30 minutes.
12. Dilute secondary antibody (usually 1:250 for lymph glands) in blocking solution at desired concentration.
13. Discard blocking solution from lymph glands and replace with 100  $\mu$ L or greater volume of primary Incubate in secondary antibody (diluted usually 1:250) in block for 3 hours at room temperature or overnight at 4C.
14. Remove secondary antibody.
15. Add 1mL 1XPBST, place on mixer for 15 minutes, discard buffer; repeat once.
16. Add 1mL 1XPBST + DNA stain (DAPI, etc.), place on mixer for 15 minutes, discard buffer.
17. Add 1mL 1XPBS, place on mixer for 15 minutes, discard buffer; repeat once (these wash steps remove detergent and excess DNA stain).

- 18.** Use pipette and dissection microscope to remove as much buffer as possible.
- 19.** Add a drop of VECTASHIELD to the lymph glands in the 1.5 mL microfuge tube.
- 20.** Store in refrigerator or begin mounting procedure.

#### Box 4 | Mounting fixed lymph glands for microscopy

Dissection of lymph glands itself can be difficult even for people with lots of practice with other tissues such as imaginal discs, however mounting them successfully after processing is just as difficult. Fewer things are more frustrating than dissecting high quality lymph gland samples, only to destroy them during the mounting process. As with dissection, many approaches will work depending on the skill level of the investigator. Below is the primary method taught to individuals in our laboratory.

1. Transfer lymph gland complexes (in VECTASHIELD) from the microfuge tube to a glass slide using a small transfer pipette or a pipette with a 200  $\mu$ L tip cut by a razor blade to enlarge the opening.
2. Place a drop of VECTASHIELD on a new slide.
3. Using forceps, move lymph gland complexes from the first slide to the second, taking care to minimize the transfer of mounting medium.
4. Once transferred, use forceps to form the VECTASHIELD into a square or rectangle shape, leaving 5-8 mm of dry area before the side edges.
5. Carefully separate the lymph gland complexes and distribute them in the rectangle of medium. Use of single tines of the forceps is recommended because surface tension of the medium will cause it to flow readily between both tines of a single forceps, which will disrupt your tissues on the slide.
6. Place one tine of one forceps underneath the posterior end of the dorsal vessel, lift slightly, and gently pull the lymph glands (and anterior structures) toward the straight, perpendicular edge of the medium where it meets the glass slide. The brain and discs create drag that straightens out the lymph glands (dorsal vessel) as you move through the medium.
7. Draw the dorsal vessel/lymph glands out past the perpendicular edge of the medium and onto the dry part of the slide. As this occurs, the volume of medium around the lymph gland decreases, causing the lymph gland and other structures to “sit down” onto the glass surface, instead of just floating around. It is critical that the tissues make contact with the slide, otherwise the samples will dislodge and become contorted upon coverslipping.
8. Repeat this process for each lymph gland around the perimeter of the medium, such that each one is well spaced and straight.
9. Separate the brain and imaginal discs from the lymph glands of each complex. This is most easily achieved by using a scissoring motion with single tines of two different forceps (avoid placing both tines of a pair of forceps into the medium). Place the tines in a “Figure X” over the dorsal vessel anterior to the primary lobes, usually at or near the location of the ring gland. Hold one tine still while drawing the other tine across the glass surface; as the tines pass each other, the dorsal vessel is severed. Push the brain complex away from the lymph gland, toward the center of the medium pool.

- 10.** Once the lymph glands and brain complexes have been separated, discard or move most of the brain complexes (and any other tissue) to a new slide. Place one brain near each corner of the medium rectangle and one or two in the center of the medium; these will serve as spacers to protect the lymph glands from severe compression by the cover slip.
- 11.** Place one short edge of an appropriately sized cover slip down onto the glass slide near, but not in, the medium. Take care to center it longitudinally. Place the tines of a forceps under the opposite end and slowly lower the cover slip onto the medium. To avoid air bubbles and misalignment, do NOT drop it onto the medium. Wait for the medium to distribute between the cover slip and the slide. Slight pressure on the cover slip with the forceps can help distribute the medium, however be careful not to compress the tissues. If more mounting medium is needed, it can be added using a pipette tip at the edge of the cover slip. Excess medium should be removed at the edge with a pipette or wicked away using a Kimwipe.
- 12.** Once the cover slip is down, seal the cover slip around the edges using nail polish. Allow to air dry for a few minutes, then store slides (horizontal is best) at 4 °C until they can be imaged.

**Box 5 | Staining for Reactive Oxygen Species (ROS) in the lymph glands***Special Materials*


- 1X Schneider's Drosophila Medium + L-Glutamine (GIBCO, cat. no. 11720)
  - Anhydrous Dimethyl Sulfoxide (DMSO) > 99.9% (Sigma-Aldrich, cat. no. 276855)
  - Dihydroethidium – special packaging (Molecular Probes, cat no. D11347)
  - VECTASHIELD Mounting Medium (Vector laboratories, cat. no. H-1000)
  - Micro Spot Plate – 3 well (Electron Microscopy Sciences, cat. no. 71561- 01)
  - 12-well microscope slides with hydrophobic barrier (Erie Scientific).
1. Dissect lymph glands in room temperature Schneider's medium. Do not use cold 1XPBS, which may inhibit respiration, thereby interfering ROS production. Because these samples are not kept on ice or immediately fixed, dissection times less than 30 minutes are optimal for limiting tissue degradation.
  2. Quickly reconstitute DHE in anhydrous DMSO (see reagent set up). Reconstituted dye solution should appear slightly pink in color; a more intense color such as purple may be indicative of oxidation of the dye.
  3. Make DHE staining solution by adding 1  $\mu$ L of the reconstituted DHE/DMSO to 1mL of Schneider's medium (in a microfuge tube) to give a final concentration of approximately 30 $\mu$ M. Vortex sample 15-30 seconds, but not more.
  4. Remove Schneider's medium from lymph glands and replace with the DHE staining solution. Incubate 5 minutes in the dark, on an orbital shaker at room temperature.
  5. Wash three times, 5-minutes each, with Schneider's medium in the dark, on an orbital shaker at room temperature.
  6. Lightly fix for 5 minutes in 7% formaldehyde/1XPBS.
  7. Wash once in 1XPBS.
  8. Immediately mount lymph glands on a glass slide in VECTASHIELD medium.
  9. Image immediately using fluorescence microscopy.

**Box 6 | Staining for Reactive Oxygen Species (ROS) in circulating cells**

1. Bleed larvae into 20  $\mu$ L Schneider's medium. Mix thoroughly and take care to remove as little liquid volume as possible when removing the carcass.
2. Let hemocytes settle 20-30 minutes in a humidified chamber at room temperature.
3. Toward the end of the settling period, reconstitute DHE and prepare DHE staining solution as described above for lymph glands.
4. Wash settled hemocytes by gently removing Schneider's medium from the edge of each well using a pipette, and replace with 20  $\mu$ L Schneider's medium.
5. Remove Schneider's medium from each well and replace with 20  $\mu$ L DHE staining solution. Incubate 5 minutes in a dark, humidified chamber.
6. Remove DHE staining solution and wash each well of cells twice, 5 minutes each, with 20  $\mu$ L Schneiders medium.
7. Lightly fix cells with 20  $\mu$ L of 4% formaldehyde/1XPBS for 5 minutes.
8. Remove fixative and wash once with 1XPBS.
9. Remove 1XPBS and add a small volume (~2  $\mu$ L) of VECTASHIELD to each well. Coverslip the slide.
10. Image immediately using fluorescence microscopy.

**Box 7 | Collection of adult hemocytes by perfusion**

1. Anesthetize adults of interest on a standard carbon dioxide pad under the dissecting microscope
2. Using forceps and microscissors, cut a small slit or hole in the posterior abdomen; this should be slightly lateral so that internal organs, genitalia, and the gut are not disrupted.
3. Place a 10  $\mu$ L drop of 1XPBS onto a well of a 12-well hydrophobic slide.
4. Using a mouth pipette, draw 1XPBS into a drawn glass capillary needle.
5. Grasp the adult fly by the wing using forceps and gently push the tip of the needle into the lateral thorax.
6. Place cut posterior end of the adult fly in or near the drop of 1XPBS and gently perfuse buffer through the adult.
7. Monitor the drop of buffer on the slide for hemolymph streaming into it; when the drop approximately doubles in size, halt perfusion and remove the adult carcass from the drop.
8. Place slide in a humidified chamber for 30 minutes.
9. Carefully remove buffer using a pipette and replace with 20  $\mu$ L of freshly prepared fixative.
10. Return to humidified chamber for 10 minutes.
11. Carefully remove fixative and replace with 20  $\mu$ L of 1XPBS, or continue with processing as described previously.



**Figure 1. A new triple-fluorescence reporter line identifies mature, progenitor, and niche cells within the hematopoietic lymph gland**

**A)** Differentiated blood cells in the cortical zone of a third-instar larval lymph gland are marked by DsRed protein (red) expression driven by the *Hemolectin* enhancer (*Hml-DsRed*). **B)** Blood progenitor cells in the medullary zone express EBFP2 protein (cyan) under the control of the domeless MESO enhancer (*dome-MESO-EBFP2*). **C)** PSC cells express EGFP protein (green) under the control of the *hedgehog* PSC enhancer (*hh-EGFP*). **D)** A combined image (Merge) showing differential marking of blood cell populations juxtaposed in the lymph gland primary lobe. DNA (blue) is labeled with TO-PRO-3 (Molecular Probes)

in each panel. **E)** Overview of the generation of the *dome-MESO-EGFP* and *dome-MESO-EBFP2* constructs. The 2.8 kb *dome-MESO* enhancer (**1**) was amplified by PCR from *Drosophila* genomic DNA, along with an added 5' CACC nucleotide sequence for Gateway® Directional TOPO® cloning into the *pENTR/D-TOPO*® entry vector (**2**; Life Technologies). Upon directional cloning of the *dome-MESO* enhancer into the Gateway entry vector (**3**), the enhancer was recombined with the *Ganesh-G2* Gateway destination vector (**4**), thereby creating the *dome-MESO-EGFP Drosophila* transformation vector (**5**). The DNA coding sequence for nuclear-localized Enhanced Blue Fluorescent Protein 2 (EBFP2-NLS) was PCR amplified from the expression vector *pLV-EBFP2-nuc* (Addgene) as an AgeI/MfeI fragment (**6**). The *dome-MESO-EBFP2* expression vector (**7**) was subsequently created by removing the *EGFP-NLS* DNA sequence from the *dome-MESO-EGFP* vector by digesting with AgeI and MfeI restriction endonucleases and then ligating in the *EBFP2-NLS* DNA sequence into the vector using the same restriction sites.

**Table 1**

Established genetic markers for the embryonic hematopoietic system

| Hematopoietic Cells         | Marker          | Genetic Reporters                                  | Antibody * | Cell Type    | References                                                                                   |
|-----------------------------|-----------------|----------------------------------------------------|------------|--------------|----------------------------------------------------------------------------------------------|
| embryonic circulating cells | <i>srp</i>      | <i>srp-gal4 srp-HEMO-gal4 srp<sup>D</sup>-gal4</i> | P          | Pan          | Milchanowski et al., 2004 [62]<br>Bruckner et al., 2004 [108]<br>Crozatier et al., 2004 [63] |
|                             | <i>gcm</i>      | <i>gcm-lacZ gcm-gal4</i>                           | P          | PL           | Bernardoni et al., 1997 [109]<br>Cho et al., 2002 [110]                                      |
|                             | <i>Iz</i>       | <i>Iz-gal4</i>                                     | M          | CC           | Lebestky et al., 2000 [52]                                                                   |
|                             | <i>crq</i>      | <i>crq-gal4</i>                                    | P          | PL           | Franc et al., 1996 [25];<br>Olofsson and Page, 2005 [111]                                    |
|                             | <i>proPO-A1</i> |                                                    | P          | CC           | Rizki et al., 1985 [40];<br>Waltzer et al., 2002 [112]                                       |
|                             | <i>CgC25</i>    | <i>Cg-gal4</i>                                     |            | PL           | Milchanowski et al., 2004 [62]; our observation                                              |
|                             | <i>Pxn</i>      | <i>Pxn-gal4</i>                                    | P          | PL           | Nelson et al., 1994 [94];<br>Stramer et al., 2005 [68]                                       |
| embryonic lymph gland       | <i>srp</i>      |                                                    | P          | Pan          | Jung et al., 2005 [54]                                                                       |
|                             | <i>odd</i>      | <i>odd-gal4 odd-lacZ</i>                           | P          | Pan          | Ward and Skeath, 2000 [113]; Jung et al., 2005 [54]                                          |
|                             | <i>Hand</i>     | <i>Hand-gal4 Hand-GFP Hand-lacZ</i>                |            | Pan          | Han et al., 2005 [114];<br>Evans et al., 2009 [115]                                          |
|                             | <i>Antp</i>     | <i>Antp-gal4</i>                                   | M          | PSC          | Mandal et al., 2007 [64]                                                                     |
|                             | <i>hth</i>      |                                                    | P          | Pan, non-PSC | Mandal et al., 2007 [64]                                                                     |
|                             | <i>col</i>      | <i>col-gal4</i>                                    |            | Pan, PSC     | Crozatier et al., 2004 [63];<br>Krzemien et al., 2007 [78]                                   |
|                             | <i>Dot</i>      | <i>Dot-gal4</i>                                    |            | Pan, PSC     | Kimbrell et al., 2002 [116];<br>Jung et al., 2005 [54]                                       |

\* P = polyclonal; M = monoclonal

**Table 2**

Established genetic markers for the larval lymph gland

| Hematopoietic Cells              | Marker                   | Genetic Reporters                                                                                                            | Antibody * | Mature cell type | References                                                                                    |
|----------------------------------|--------------------------|------------------------------------------------------------------------------------------------------------------------------|------------|------------------|-----------------------------------------------------------------------------------------------|
| Pan                              | <i>srp</i>               |                                                                                                                              | P          |                  | Lebestky et al., 2000 [52]                                                                    |
|                                  | <i>odd</i>               | <i>odd-lacZ</i>                                                                                                              | P          |                  | Jung et al., 2005 [54]                                                                        |
|                                  | <i>He</i>                | <i>He-gal4</i>                                                                                                               | M          |                  | Kurucz et al., 2003 [117]                                                                     |
| Posterior Signaling Center (PSC) | <i>Antp</i>              | <i>Antp-gal4</i>                                                                                                             | M          |                  | Mandal et al., 2007 [64]                                                                      |
|                                  | <i>Ser</i>               | <i>Ser9.5-lacZ</i><br><i>Ser9.6-gal4</i>                                                                                     | P          |                  | Lebestky et al., 2000 [52]<br>Jung et al., 2005 [54]                                          |
|                                  | <i>hh</i>                | <i>hh-GFP</i>                                                                                                                | P          |                  | Mandal et al., 2007 [64]; Tokusumi et al., 2009 [81]                                          |
|                                  | <i>col</i>               | <i>col-gal4</i>                                                                                                              |            |                  | Crozatier et al., 2004 [63]; Krzemien et al., 2007 [78]                                       |
| Progenitors/Medullary Zone       | <i>Pvf1</i>              |                                                                                                                              | P          |                  | Mondal et al., 2011 [84]                                                                      |
|                                  | <i>dome</i>              | <i>dome-gal4</i><br><i>dome-gal4; ELAV-gal80</i><br><i>dome-MESO-lacZ</i><br><i>dome-MESO-EGFP</i><br><i>dome-MESO-EBFP2</i> |            |                  | Jung et al., 2005 [54]<br>Krzemien et al., 2007 [78]<br>Mondal et al., 2011 [84]<br>this work |
|                                  | <i>DE-Cad</i>            |                                                                                                                              | M          |                  | Jung et al., 2005 [54]                                                                        |
|                                  | <i>upd3</i>              | <i>upd3-gal4</i>                                                                                                             |            |                  | Jung et al., 2005 [54]                                                                        |
|                                  | <i>ptc</i>               |                                                                                                                              | M          |                  | Mandal et al., 2007 [64]                                                                      |
|                                  | <i>wg</i>                |                                                                                                                              | M          |                  | Sinenko et al., 2009 [85]                                                                     |
|                                  | <i>ci</i>                |                                                                                                                              | M          |                  | Mandal et al., 2007 [64]                                                                      |
|                                  | ROS <sup>high</sup>      |                                                                                                                              |            |                  | Owusu-Ansah and Banerjee, 2009 [89]                                                           |
|                                  | <i>bam</i>               |                                                                                                                              | M          |                  | Tokusumi et al., 2011 [88]                                                                    |
|                                  | p-CamK-II                |                                                                                                                              | P          |                  | Mondal et al., 2011 [84]                                                                      |
|                                  | p-Akt                    |                                                                                                                              | P          |                  | Shim et al., 2012 [87]                                                                        |
|                                  |                          |                                                                                                                              | P          |                  | Shim et al., 2013 [118]                                                                       |
|                                  | GABA-R <sup>high</sup>   |                                                                                                                              | P          |                  | Shim et al., 2013 [118]                                                                       |
|                                  | <i>TepIV</i>             | <i>TepIV-gal4</i><br><i>TepIV-gal4; ELAV-gal80</i>                                                                           |            |                  | Irving et al., 2005 [101]; Krzemien et al., 2007 [78]<br>Avet-Rochex et al., 2010 [86]        |
|                                  | <i>Pvr<sup>low</sup></i> |                                                                                                                              | P          |                  | Jung et al., 2005 [54]; Mondal et al., 2011 [84]                                              |

| Hematopoietic Cells                 | Marker                      | Genetic Reporters                                  | Antibody * | Mature cell type | References                                                                   |
|-------------------------------------|-----------------------------|----------------------------------------------------|------------|------------------|------------------------------------------------------------------------------|
| Differentiating cells/Cortical Zone | <i>NimC1 (P1)</i>           |                                                    | M          | PL               | Kurucz et al., 2007 [33]; Kurucz et al., 2007 [96]                           |
|                                     | <i>Hml</i>                  | <i>Hml-gal4</i><br><i>Hml-DsRed</i>                |            | PL               | Sinenko et al., 2004 [119]; Makhijani et al., 2011 [58]                      |
|                                     | <i>eater</i>                | <i>eater-gal4</i><br><i>eater-GFP</i>              |            | PL               | Tokusumi et al., 2009 [81]                                                   |
|                                     | <i>Pvr<sup>high</sup></i>   |                                                    | P          |                  | Jung et al., 2005 [54]; Mondal et al., 2011 [84]                             |
|                                     | <i>Pxn</i>                  | <i>Pxn-gal4</i>                                    | P          | PL               | Jung et al., 2005 [54]; Stofanko et al., 2008 [60]                           |
|                                     | <i>Collagen (gal4)</i>      | <i>Cg-gal4</i>                                     | M          | PL               | Jung et al., 2005 [54]                                                       |
|                                     | <i>GABA-R<sup>low</sup></i> |                                                    | P          |                  | Shim et al., 2013 [118]                                                      |
|                                     | <i>Iz</i>                   | <i>Iz-gal4</i>                                     | M          | CC               | Lebestky et al., 2000 [52]; Jung et al., 2005 [54]                           |
|                                     | <i>proPO-A1</i>             | <i>Bc<sup>1</sup></i><br><i>Bc-GFP, -BFP, -RFP</i> | P          | CC               | Rizki et al., 1980 [120]; Jung et al., 2005 [54]; Tokusumi et al., 2009 [81] |
|                                     | <i>Sima</i>                 |                                                    | P          | CC               | Muhkerjee et al., 2011 [121]                                                 |
|                                     | <i>ItgaPS4</i>              |                                                    | P          | LM               | Crozatier et al., 2004 [63]; Irving et al., 2005 [101]                       |
|                                     | <i>msn</i>                  | <i>msn-lacZ</i><br><i>msn-GFP, msn-RFP</i>         |            | LM               | Braun et al., 1997 [98]; Tokusumi et al., 2009 [99]                          |
|                                     | <i>atilla (L1)</i>          | <i>atilla-GFP</i>                                  | M          | LM               | Honti et al., 2009 [102]                                                     |
|                                     | <i>mys</i>                  |                                                    | M          | LM               | Irving et al., 2005 [101]                                                    |

\* P = polyclonal; M = monoclonal