

# UC Irvine

## UC Irvine Previously Published Works

### Title

Purification of Human Monocytes and Lymphocyte Populations by Counter Current Elutriation- A Short Protocol.

### Permalink

<https://escholarship.org/uc/item/6hx375rz>

### Journal

Bio-protocol, 3(23)

### ISSN

2331-8325

### Authors

Clarke, Elizabeth V  
Benoit, Marie E  
Tenner, Andrea J

### Publication Date

2013

### DOI

10.21769/bioprotoc.981

Peer reviewed



Published in final edited form as:

*Bio Protoc.* 2013 December 5; 3(23): .

## Purification of Human Monocytes and Lymphocyte Populations by Counter Current Elutriation– A Short Protocol

Elizabeth V. Clarke, Marie E. Benoit, and Andrea J. Tenner\*

Molecular Biology and Biochemistry Department, University of California, Irvine, CA, USA

### Abstract

Investigations of the activation processes involved in human monocytes and monocyte-derived macrophages and dendritic cells often required large numbers of cells that have not been possibly altered or activated by adherence to surfaces, by binding of antibodies to surface antigens during positive selection, or by release of activators by platelets or other non myeloid cells during isolation or co-culture. Human peripheral blood monocytes as well as lymphocytes from the same blood donor can be isolated by counterflow elutriation using a modification of the technique of Lionetti *et al.*, 1980 as described previously (Bobak *et al.*, 1986). From a unit of blood drawn into anticoagulant, 60-120 million monocytes can be obtained. These cells are not activated and have been shown to be appropriately capable of differential activation in multiple studies.

### Materials and Reagents

1. Peripheral blood from normal human donor (450 ml)
2. CPD AS-1 blood collection bags (Citrate-Phosphate-Dextrose, Adenine saline-1) (Baxter)
3. 1× PBS-2mM EDTA
4. Sterile milliQ H<sub>2</sub>O
5. Lymphocyte separating medium (LSM) (MP Biomedicals, catalog number: 0850494)
6. HBSS (Fisher Scientific, catalog number: 21023CV)
7. Ethanol (70%, 1 liter)
8. 10% Bleach (1 liter)
9. 10% defined FBS (Hyclone, catalog number: SH30070.03)
10. 1% P/S (Mediatech, CellGro<sup>®</sup>, catalog number: 30-002-CIRF)
11. 1% L-Gln (Mediatech, CellGro<sup>®</sup>, catalog number: 25-005-CI)
12. Recombinant human IL-2 (Pepro Tech, catalog number: 200-02)
13. Recombinant human M-CSF (Pepro Tech, catalog number: 300-25)

\*For correspondence: atenner@uci.edu.

14. 25% Albumin (Human), USP (Talecris Plasma Resources, catalog number: 13533-684-20)
15. Complete media (see Recipes)
16. 1× sterile PBS (dilute 1:5 from 5× PBS pH 7.0 – 7.2) (see Recipes)
17. Ammonium Chloride Potassium (ACK) buffer (see Recipes)
18. 1× sterile Elutriation Buffer (see Recipes)
19. 5× Phosphate Buffered Saline (see Recipes)

## Equipment

1. Extruder
2. T75 flasks
3. Nunc blue cap 50-ml polypropylene tubes (important as monocytes stick to other tubes lowering yields)
4. Tubing clamps
5. Scissors
6. P1000 Pipetteman
7. Petri dishes
8. 250 ml plasma tube
9. Isopropanol wipes, individually packaged
10. Beckman JE6B Elutriation System and Rotor (Beckman Coulter, model: JE6B-IM-3) (November 1992)\*\*
11. Centrifuge (Sorvall, model: RT6000 or RC-3B)
12. Stericup (Fisher Scientific, catalog number: SCGVU11RE)

## Procedure

### A. Receive 450 ml peripheral blood from normal human donor drawn in CPD AS-1

Using an RC-3B Centrifuge, spin a 450 ml bag of blood upright with no brake, 25 °C, 1665 × *g*, for 7 min. Takes approximately 30 min to come to a stop without the brake.

### B. Assemble elutriator and chambers

See Beckman JE6B Elutriation System and Rotor Assembly Manual.

### C. Sterilize and prepare elutriation chamber

1. Attach 1 ml pipet to elutriation pump tubing and transfer tubing (by submerging in each solution) to each of the following in sequence:

---

\*\*Beckman Coulter has newer models Avanti J-301 or J-26S XP with Elutriation Systems, but we have no experience with those.

- a. Wash #1 10% fresh bleach (in sterile milliQ H<sub>2</sub>O)
  - b. Wash #2 300-500 ml: sterile milliQ H<sub>2</sub>O
  - c. Wash #3 300-500 ml: 70% EtOH (sterile milliQ H<sub>2</sub>O)
  - d. Spin elutriator rotor by hand for 30 sec to remove air pockets.
  - e. Wash #4 300-500 ml: sterile milliQ H<sub>2</sub>O
  - f. Wash #5 300-500 ml: sterile PBS
2. Finally, equilibrate tubing and chamber with elutriation buffer (with 100-200 ml approximately 11-12 ml/min.flow rate).

#### D. Collect blood from centrifuge

1. Sit blood bag on extruder carefully and place onto the hooks on extruder.
2. Wipe tube of blood bag and scissors with alcohol wipes to sterilize. Wipe tubing end with alcohol, attach clamp, and insert into a 250 ml plasma tube (if collecting plasma for later use) or waste (if not collecting plasma).
3. With clamp shut on tubing, release press to squeeze bag gently.
4. Open clamp to control flow to a fast drip/slow stream. Stop plasma collection when plasma is approximately 1 inch from top of bag.
5. Transfer blood bag tube to T75 flask containing 50 ml PBS 0.02 mM EDTA, where you will collect approximately 25 ml plasma and 20 ml erythrocytes plus white cells) often called the buffy coat. Swirl flask gently. Total volume in flask should be approximately 90-100 ml.
6. Tilt 50 ml tube containing 14-18 ml LSM and using a 10 ml pipette, layer approximately 30 ml buffy coat-PBS-EDTA carefully from flask onto LSM, repeat with second and third LSM tube.
7. Remainder of buffy coat-PBS-EDTA should be layered onto last LSM tube.
8. Balance LSM tubes, spin in RT6000 at  $350 \times g$  for 40 min, 22 °C with brake off.

#### E. Blood/cells

1. Remove diluted Plasma from top of tubes of LSM (approximately 25 ml).
2. Remove peripheral blood mononuclear cells (PMBC) interfaces, add to 40 ml PBS/1 mM EDTA in 50 ml tubes and invert to mix.
3. Spin in RT6000 to remove platelets,  $200 \times g$  for 10 min. Pellet should be visible.
4. Pour off most of the supernatant.
5. Resuspend all pellets, pooling in a total of 8–10 ml elutriation buffer. Avoid air bubbles.
6. Count PBMCs: should be between  $7 \times 10^8$  and  $1.5 \times 10^9$  total cells depending on the donor.

## F. Elutriation

1. Run some elutriation buffer through to equilibrate (100 ml).
2. Spin rotor manually for approximately 30 sec to remove air pockets.
3. Close lid and allow speed to get up to  $1450 \times g$  (final speed) before checking flow rate: should be between 11.4 and 11.6 ml/min.
4. Transfer pipette attached to pump tubing into cells and pump full volume of cells (8-10 ml) before transferring pipette back to elutriation buffer. Avoid taking up any air by pinching tubing during transfer.
5. Collect the effluent. Collect the first two 50 ml tubes of effluent which will contain the lymphocyte population (with some contaminating platelets and red cells) and spin down at  $100 \times g$  to remove platelets. Lyse any remaining erythrocytes in ACK at room temperature (by resuspending the pellet in the ACK and gently inverting the tube 1-2 times), wash the lymphocytes once in HBSS by resuspending in 25 ml, count cell concentration, pellet cells again at  $100 \times g$ , and resuspend the lymphocytes in complete media containing 50 U/ml IL-2 at 2 million/ml in complete media in T75 flasks. Monocytes will remain in chamber, due to size and density and centrifugation conditions.
6. After 35 min run at 11.5 ml/min increase speed to 12.5 ml per min. Continue to run for 25 min (total runtime is 60 min).

## G. End of run

1. Clamp output and input tubes before pump is turned off, then immediately stop centrifuge.
2. Open lid and unlock black lock, then remove input tube first, then lower outlet tube.
3. Remove elute chamber, with spout up.
4. In hood, using a P1000 remove the 6.3 ml of sample from the elutriation chamber and add to 50 ml tube.
5. Bring up volume up to 20 ml with  $1 \times$  PBS. Count cells. Monocytes should be approximately 10-20% of original PBMC count. Wash cell pellet with sterile PBS.  $200 \times g$  for 10 min. Resuspend pellet to  $100 \times 10^6$  cells/ml.
6. Plate monocytes in 10 cm petri dishes at  $0.5 \times 10^6$ /ml (10 ml total/dish) in complete media containing 25 ng/ml M-CSF for derivation into macrophages.

## H. Clean elutriator

1. Place elute chamber back into rotor and rotor back into centrifuge chamber attaching tubes to wash out/decontaminate rotor, tubing and chamber.
2. Undo clamps and run through 300 ml 10% bleach, 300 ml autoclaved H<sub>2</sub>O, 300 ml 70% EtOH, 300 ml water again. Disassemble rotor, set in 70% ethanol for 30% and let air dry overnight.

## Recipes


1. Complete media
  - RPMI
  - 10% defined FBS
  - 1% P/S
  - 1% L-Gln
2. ACK buffer
  - To 450 ml milliQ water, add:
  - 4.145 g Ammonium Chloride (annhydrous)
  - 0.5 g potassium bicarbonate
  - 18.6 mg disodium EDTA
  - Adjust pH to 7.4
  - Bring final volume to 500 ml with milliQ water
  - Filter sterilize (Using Stericup) and store at 4 °C.
3. 1× sterile Elutriation Buffer, pH 7.3
  - 8.3 g NaCl
  - 5 g dextrose in 1 L MQ H<sub>2</sub>O
  - 10 mM Na<sub>2</sub>HPO<sub>4</sub>
  - 3 mM NaH<sub>2</sub>PO<sub>4</sub>\*H<sub>2</sub>O
  - 100,000 units Pen/Strep
  - 0.625% Human Serum Albumin (HSA) from Talecris Therapeutics
4. 5× Phosphate Buffered Saline (PBS for washing cells), pH 7.0-7.2
  - 9.51 g dibasic phosphate, anhydrous (diluted is 6.7 mM Na<sub>2</sub>HPO<sub>4</sub>)
  - 4.56 g monobasic phosphate (diluted is 3.3 mM Na H<sub>2</sub>PO<sub>4</sub>-H<sub>2</sub>O)
  - 85 g NaCl (diluted will be .14 M)
  - Dissolve and bring to final volume of 2000 ml with MilliQ (LPS-free) water.
  - Filter sterilize.
  - Diluted at room temperature should be pH 7.0-7.2, and K = 14-16. Store at 4 °C.

## Acknowledgments

The original protocol was published in Bobak *et al.* (1986). Further modification and development of these methods was supported in part by NIH AI-41090 and T32 AI 60573.

## References

1. Bobak DA, Frank MM, Tenner AJ. Characterization of C1q receptor expression on human phagocytic cells: effects of PDBu and fMLP. *J Immunol.* 1986; 136(12):4604–4610. [PubMed: 3486907]
2. Lionetti, F.; Hunt, S.; Valeri, C. *Methods of Cell Separation.* Plenum Publishing Corporation; New York: 1980.



**Figure 1.**

Adapted from <http://www.genex.cl/stock/Mediatech.html>